

BEDSIDE TABLE - height adjustable tray - wood for 27672/3



Product Code	27676
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1391274
CND	V0899
EAN/UPC	8023279276763
UMDNS	10-531

Description

WOODEN BEDSIDE TABLE - height adjustable, tiltable tray

Manufactured by quality laminated wooden panel.
Particularly functional: with drawer, 2 doors (can be used from both sides) and lateral height adjustable (73-100 cm) withj tiltable tray at 90°.

Size: 50 x 33 x h 85 cm.
Lateral tray: 58 x 33 cm.
4 revolving castors (2 lockable).
ABS border.