

# N-SHAPED BEDRAIL 165.5x40 cm



|               |                |
|---------------|----------------|
| Product Code  | 27697          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 2518291       |
| CND        | V08068001     |
| EAN/UPC    | 8023279276848 |

## Description

N-SHAPED BEDRAIL 165x40 cm

Bedrail made of steel frame and surface with coating.  
Siderail can be fixed by screw.  
Max load: 75 kg.

Multilingual manual: GB, FR, IT, ES, PT, DE, PL, RO, SE, GR.