

COMMODE CHAIR - stainless steel



Product Code	27701
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2184712
CND	Y091203
EAN/UPC	8023279277012
PARAF	903118255

Description

COMMODE CHAIR - chromed steel

Chromed steel chair frame, with plastic detachable commode seat.

Technical Specifications

- Total width: 56 cm
- Seat width: 44 cm
- Seat height: 47 cm
- Total height: 94 cm
- Seat depth: 43 cm
- Seat height: 38 cm
- Max. load: 100 kg
- Kit size: 50 x 19 x 75 cm
- GW (kg): 8.5