



Product Code	27703
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2184713
CND	Y122199
EAN/UPC	8023279277036
PARAF	903607366

Description

COMMODE WHEELCHAIR - stainless steel

Stainless steel chair frame, detachable armrest and footrest, rear castors with lock, plastic commode seat.

Multilingual manual: GB, FR, IT, ES.

Technical Specifications

- Total width: 48 cm
- Seat width: 44 cm
- Rear wheel diameter: 12.5 cm
- Front wheel diameter: 12.5 cm
- Seat height: 51 cm
- Total wheelchair height: 92 cm
- Total length: 90 cm
- Seat depth: 44 cm
- Seat height: 42 cm
- Max. load: 100 kg
- Kit size: 54 x 53 x 53 cm
- NW GW (kg): 17-18