

STANDARD FOLDING WHEELCHAIR



Product Code	27709
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2184716
CND	Y122106
EAN/UPC	8023279277098
UMDNS	16949
PARAF	903607683

Description

STANDARD WHEELCHAIR - 46 cm seat - blue tissue

- single cross bar
- chromed frame
- fixed armrest, PVC pad
- fixed footrest, flip-up plate with leg strap
- 20 cm (8") solid front castors
- 61 cm (24") solid rear wheels
- locking brakes on rear wheels
- large pocket on the back

Technical Specifications

- Total width: 64 cm
- Folded wheelchair width: 24 cm
- Seat width: 46 cm
- Rear wheel diameter: 58 cm
- Front wheel diameter: 20 cm
- Seat height: 49 cm
- Total wheelchair height: 87 cm
- Total length: 106 cm
- Seat depth: 40 cm
- Seat height: 45 cm
- Max. load: 100 kg
- Kit size: 94 x 23 x 89 cm
- Net Weight: 17 Kg
- Gross Weight: 19.2 Kg