

EXTRALARGE WHEELCHAIR - Steel



Product Code	27716
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2184716
CND	Y122106
EAN/UPC	8023279277166
UMDNS	16949
PARAF	977812015

Description

EXTRA LARGE WHEELCHAIR - 56 cm seat - black tissue

- double cross bar
- powder coated steel frame
- super width 56 cm seat
- flip-up armrest, PU pad
- detachable footrest, flip-up plate with leg strap
- 20 cm (8") rubber foam front castors
- 61 cm (24") pneumatic rear wheels
- large pocket on the back

Technical Specifications

- Total width: 73 cm
- Folded wheelchair width: 32 cm
- Seat width: 56 cm
- Rear wheel diameter: 62 cm
- Front wheel diameter: 20 cm
- Seat height: 54 cm
- Total wheelchair height: 93 cm
- Total length: 109 cm
- Seat depth: 45 cm
- Seat height: 46 cm
- Max load: 125 kg
- Kit size: 83 x 33 x 93 cm
- Net Weight: 21.5 Kg
- Gross Weight: 24.1 Kg