

DELUXE WHEELCHAIR - Aluminium



Product Code	27717
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2184715
CND	Y122106
EAN/UPC	8023279277173
UMDNS	16949
PARAF	977812027

Description

DELUXE WHEELCHAIR - 46 cm seat - black tissue

- double cross bar
- flip-up height adjustable armrest, PU pad
- detachable footrest, flip-up plate with leg strap
- 61 cm (24") solid PU rear wheels
- 15 cm (6") solid front castors
- locking brakes on rear wheels
- anti tipper system with 2 small wheels

Large pocket on back
Removable padding with velcro fixing
Leg strap

Technical Specifications

- Total width: 67 cm
- Folded dimensions: 36 cm
- Seat width: 46 cm
- Rear wheel diameter: 61 cm
- Front wheel diameter: 15 cm
- Seat height: 51 cm
- Total wheelchair height: 91 cm
- Total wheelchair length: 105 cm
- Backrest height: 40 cm
- Seat depth: 45 cm
- Max load: 100 kg
- Kit size: 83 x 28 x 93 cm
- Net Weight: 17 Kg
- Gross Weight: 19 Kg