

# LEONARDO WOOD STICK - metal engraved "T" handle



Product Code	27736
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8050714360541

## Description

BEECHWOOD STICKS:

LEONARDO METAL ENGRAVED "T" HANDLE - unisex - black

Max load for wooden sticks: 100 kg.

Height: 92 cm (Can change from 90 to 96 cm depending on rubber tip and handle)

Made in Italy.