

TRUBOND NON ABSORB. SUTURE gauge 5/0 circle 3/8
needle 20 mm - 76 cm - green



Product Code	22774
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2380199
CND	H010201010202
EAN/UPC	8023279227741

Description

TRUBOND - non absorbable surgical suture, silicone coated and braided polyester for permanent support and strength

Equivalent Ethicon brand: ETHIBOND/MERSILENE

- Gauge: 5/0
- Metric: 1
- Colour: Green
- Circle: 3/8
- Needle profile: Round body double arm
- Needle length: 20 mm
- Suture length: 76 cm

- Technical details:
- braided polyester, coated with silicone for a smooth knot tie down
 - high tensile strength for reliability
 - low tissue reactivity - undergoes no significant change in body
 - excellent handling and long lasting
 - sterilized by ethylene oxide
 - shelf life of 5 years

- Applications:
- used for permanent and removable fixtures
 - prosthetic implantations
 - "Face lifts"
 - cardiovascular surgery
 - ideal for valve repair/ replacement

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.