

CONTACT PLUS SOFT RAISED TOILET SEAT with lid



Product Code	27745
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2044555
CND	Y091203
EAN/UPC	3701127704433
PARAF	987471721

Description

CONCTACT PLUS SOFT RAISED TOILET SEAT - with lid

High quality soft toilet raiser, made of high-density foam injected in a waterproof PVC skin.

It fits most toilet bowl in a very secure way.

This comfortable seat with ergonomic shape and high density foam, raises the toilet by 11 cm (4.3 inches).

Tested for 185 kg (407 lb).

Dimensions outside: 37 x 41 x h 16 cm

Dimensions inside: 21 x 26 x h 16 cm

Made in France.