

RAISED TOILET SEAT with fixed armrest - height 12,5 cm



Product Code	27750
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2566853
CND	Y091215
EAN/UPC	8023279277500
PARAF	903118281
GTIN	08023279277500
GMDN	31059

Description

RAISED TOILET SEAT height 12.5 cm

With fixing system (front screw) and fixed armrest made in anodized aluminium.
Max load 130 kg (286.6 lbs)
Packing: 410 x 165 x 468 mm

Fixed armrest.

Instructions: GB, IT, FR, ES, DE, GR

Made in Italy.