



Product Code	27752
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2272691
CND	V0805030102
EAN/UPC	8023279277524
PARAF	987472659

## Description

Lifting movement through a hydraulic pump with electrical controls.

Up and down movements through an electric motor 24V D.C. and push button.

Battery fed with current transformer with an emergency push-button for the immediate stop of the movement and free up/down piston-stroke.

Painted steel frame. Anti trauma push-handle. 2 rear castors Ø 100 mm, with independent brake and 2 front castors Ø 70 mm.

Supplied with universal sling (27767).

Made in Italy.