

LEONARDO WOOD STICK - metal - curved handle



Product Code	27768
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279277685
PARAF	903607998

Description

BEECHWOOD STICKS:

LEONARDO METAL CURVED HANDLE - unisex - black

Max load for wooden sticks: 100 kg.
Height: 92 cm (Can change from 90 to 96 cm depending on rubber tip and handle)

Made in Italy.