

PATIENT TROLLEY - REMOVABLE TOP



Product Code	27801
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2831159
CND	V0899
EAN/UPC	8023279278019
PARAF	903113468
GTIN	08023279278019
GMDN	35892

Description

PATIENT TROLLEY - removable top

Made of chrome plated steel tube. Comes complete with mattress made of washable self extinguishing plastic material.
Mounted on 4 castors Ø 100 mm two of which with brakes.
Size: 180x54xh 80 cm.
Load: 110 kg.

Adjustable backrest up to a 60°

Multilingual instructions: GB, FR, IT, ES

Made in Italy.