

I.V.STAND ON 5 WHEELS TROLLEY with handle and shelf - stainless steel - 4 hooks



Product Code	27815
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2879920
CND	V08068099
EAN/UPC	8023279278156

Description

Stainless steel 201 IV stand

Adjustable height: 110-200 cm

Weight: 3.6 kg

Base material: Stainless steel

Hook (number): 4

Hook material: Stainless steel

Load kg Per hook: 2 kg

Total load: 8 kg

Castors with brake: 5

Castors Antistatic: -