

I.V. INFUSION STAND ON 5 WHEELS TROLLEY with handle and shelf - s/s - 4 hooks



Product Code	27816
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Not required
CH-REP	Not required

RDM (NSIS)	2879920
CND	V08068099
EAN/UPC	8023279278163

Description

V STAND FOR INFUSION PUMPS

IV stand made of high quality 201 stainless steel.

Lightweight and easy to move, suitable for supporting single or dual-channel infusion pumps.

Slip resistant and powder coated steel base with ABS dust cover to prevent dust accumulation, ensuring more hygienic environment. It is equipped with a multifunctional aluminium tray with integrated handle, which can support infusion pump; four stainless steel hooks and five castors Ø 76 mm, two of which with brakes.

Department name card for easy identification.

Ø 50 cm base.

4 S/S hooks

Multifunctional tray

Technical Specifications

Height adjustment: 1500 - 2500 cm

Weight: 6.4 kg

Tray size: 20x17.5 cm

Tray load capacity: 10 kg