



Product Code	27818
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	No

RDM (NSIS)	2804663
CND	Z129099
EAN/UPC	8023279278187

## Description

Stainless steel AISI 202 IV stand

Adjustable height: 135-215 cm

Weight: 3.5 kg

Base material: Aluminium

Hook (number): 4

Hook material: Stainless steel

Load kg Per hook: 3 kg

Total load: 12 kg

Castors with brake: 2

Castors Antistatic: 5