



Product Code	27829
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	470777
CND	V0899
EAN/UPC	8023279278293

Description

Infusion stand with hydraulic pump

Made of chrome plated steel.
Four chrome plated steel hooks.
Five spoke stainless steel base Ø 60 cm with castors Ø 50 mm, two with brake.
In conformity with MDR.

Made in Italy.

Technical Specifications

Adjustable height (cm): 189-260
Load (kg): 30 (5 kg each hook, 10 kg for infusion pump frame)
Weight (kg): 18
Infusion pump frame (cm): 39 x 75
Packaging size (cm): 185 x 18 x 43