



|               |                   |
|---------------|-------------------|
| Product Code  | 25783             |
| Unit of sale  | box of 50 pcs     |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8013183130213 |
| PARAF   | 977810845     |

Description

CS line sharps containers  
• CS 3 I - 16.7x16.7xh 23.2 cm

Full range of Sharps containers for the disposal of sharp or cutting waste.  
Equipped with temporary and definitive closure, certified hermetic seal, do not release any harmful emission during incineration.  
The containers hold the UN mark for the transportation of dangerous materials.

- practical opening system
- temporary and permanent closure
- needles removal
- insulin pen removal device
- maximum filling level
- shock resistant
- puncture resistant
- easy to use
- easy to assemble
- excellent stability
- square or round shaped
- sturdy handle for easier transport
- customisable label for quantities
- autoclavable at 134°C for 18'
- no emission of any harmful gas during the incineration process

Made in Italy.