

I.V.STAND ON 5 WHEELS TROLLEY - 2 hooks - disassembled



Product Code	27840
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1445594
CND	V0899
EAN/UPC	8023279278408

Description

Epoxy painted IV stand

Adjustable height: 160-220 cm
Weight: 4.0 kg
Base material: Painted steel
Hook (number): 2
Hook material: Plastic
Load kg Per hook: 5 kg
Total load: 10 kg
Castors with brake: -
Castors Antistatic: -