

# I.V.STAND ON 5 WHEELS TROLLEY - stainless steel - 4 plastic hooks



|               |                |
|---------------|----------------|
| Product Code  | 27844          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No  |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2823532        |
| CND        | Z129099        |
| EAN/UPC    | 8023279278446  |
| PARAF      | 977812179      |
| GTIN       | 08023279278446 |
| GMDN       | 36069          |

## Description

Stainless steel AISI 304 IV stand

Adjustable height: 135-215 cm

Weight: 2.9 kg

Base material: Stainless steel

Hook (number): 4

Hook material: Plastic

Load kg Per hook: 1.2 kg

Total load: 4.8 kg

Castors with brake: 1

Castors Antistatic: 5