

I.V BOTTLE HOLDER - stainless steel



| | |
|---------------|----------------|
| Product Code | 27848 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | No |

| | |
|------------|----------------|
| RDM (NSIS) | 2805613 |
| CND | Z129099 |
| EAN/UPC | 8023279278484 |
| GTIN | 08023279278484 |

Description

Stainless steel AISI 304 bottle holder for bottle > or = 1 l
Capacity 2.5 kg