



Product Code	27850
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2831159
CND	V0899
EAN/UPC	8023279278507
PARAF	903113569
GTIN	08023279278507
GMDN	35892

Description

TRENDELENBURG TROLLEY

Made of steel-alloy, oven painted ral 9010, with chrome plated handles.
Adjustable height from 50 to 70 cm, with oleodinamic pump.
Trendelenburg up to 9°.
Anti-Trendelenburg up to 14° (with gas pump).
Plastic bumpers at four corners. Maximum load: 120 kg.
4 wheels of Ø 200 mm (2 with brake).
Size: 65 x 184 (internal) x h 50-70 cm, 72 x 191 cm (external).

Steel-alloy
Ø 200 mm wheels
Oleodinamic pump for height adjustment
Bumpers at all 4 corners

Made in Italy.