



Product Code	27857
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2831153
CND	V0899
GTIN	08023279278576
GMDN	36069

Description

HYPODERMOCLYSIS STAND

Stand with hydraulic height adjustment for hypodermoclysis, for washing and surgical irrigation. Hydraulic pump with pedal for height adjustment and gripping handle for easy transport.

LOAD 28 KG

4 stainless steel AISI 304 hooks, with antirelease safety device

Stainless steel adjustable gripping handle for easy transport

Stainless steel AISI 304 stand Ø 60 mm

Pedal driven shock absorbing pump.

Height adjustable: 147-247 cm

Weight: 22 kg

Stainless steel base with 5 pivoting castors Ø 75 mm (2 antistatic with brake, 3 standard) and grey vipla bumpers.

Made in Italy.