

DIATERMO CART P with power - medium



Product Code	27866
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	509652
CND	V0880
EAN/UPC	8023279278668

Description

DIATERMO CART for electrosurgical units with power - medium

Ergonomic trolley designed to hold and carry electrosurgical units (including Argon units).
The height of 915 mm with inclined top shelf allows easy reading of the unit's display.
The trolley is provided with a second shelf with antiskid raised edges, two drawers placed at the base and a rear space with securing straps to hold two Argon cylinders.
Upright stands can hold optional accessories.
Made of epoxy-painted, tempered steel and aluminium.
With 4 socket group, fuses box, independent grounding plug, power cable 3m.

Made in Italy.

Technical Specifications

Load capacity: - Total load 200 kg - Inclined shelf 50 kg
Inclined shelf: 42x45 cm
Drawers and shelf: 42x36 cm