



Product Code	27867
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	289319
CND	V0880
EAN/UPC	8023279278675

Description

DIATERMO CART for electro-surgical units with power - large

Ergonomic trolley designed to hold and carry electro-surgical units (including Argon units).

The height of 915 mm with inclined top shelf allows easy reading of the unit's display.

The trolley is provided with a second shelf with antiskid raised edges, two drawers placed at the base and a rear space with securing straps to hold two Argon cylinders.

Upright stands can hold optional accessories.

Made of epoxy-painted, tempered steel and aluminium.

With 4 socket group, fuses box, independent grounding plug, power cable 3m.

Made in Italy.

Technical Specifications

Load capacity: - Total load 200 kg - Inclined shelf 50 kg

Inclined shelf: 50x36 cm

Drawers and shelf: 50x36 cm