

SUPER EASY CART with shelf and monitor support



Product Code	27869
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	173922
CND	V0880
EAN/UPC	8023279278699
GMDN	173922

Description

SUPER EASY CARTS - with handle

- Cart + shelf + base + monitor support (Can not be used with 27874)
- Shelf size (cm): 40x36

A new stylish and quality cart you can build according to your needs. Carts with base and 1 or 2 sturdy shelves, aluminium columns and 4 twinned castors Ø 65 mm (2 with brake). Illustrated assembly instructions.

BUILD UP YOUR CART, YOU CAN ADD: EXTRA SHELF, BASKET, INFUSION STAND, CARDIO ARM.

Made in Italy.

Technical Specifications

Dimensions: 47x42x102 cm
Load capacity: total 120 kg - shelf 10 kg
4 twinned castors Ø 65 mm (2 with brake)
Material:
- extrusions in aluminium
- plastic parts in ABS
Shelf: steel sheet epoxy powder coated