

VALUE CABINET - 1 door - tempered glass



Product Code	27898
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279278989

Description

VALUE CABINET - 1 door

Cabinet made of enameled steel sheet with tempered glass door.
Hinged glass door with lock and key.
Three stainless steel shelves adjustable in height.
Delivered in kit form.

Size: 52x35x138 cm
Pack size: 65x26x134 cm
Weight: 31 kg

Made in Italy.