

VALUE CABINET - 4 hinged doors - tempered glass



Product Code	27900
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279279009

Description

VALUE CABINET - 4 doors (hinged)

Cabinets made of enamelled steel sheet with 2 overlapped compartments.
Upper compartment with 2 tempered glass hinged doors, lock with key and 2 stainless steel, height adjustable shelves. Inferior compartment with 2 hinged doors made of enamelled steel plate and 1 enamelled steel sheet shelf.
Delivered in kit form.

Size: 75x37x164 cm
Pack size: 72x24x156 cm
Weight: 50 kg

Made in Italy.