

CABINET - 1 door - tempered glass



Product Code	27901
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	Not required

EAN/UPC	8023279279016
PARAF	903608267
GTIN	08023279279016

Description

Cabinet with epoxy coated steel structure and tempered glass.
Rack-mounted adjustable glass shelves.
Rubber feet.

Made in Italy.

Technical Specifications

- Doors: 1
- Adjustable shelves: 3
- Weight: 30 kg
- Size: 53 x 36 x h 144 cm