

CABINET - 4 doors - tempered glass



Product Code	27903
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	Not required
EAN/UPC	8023279279030
PARAF	903608293
GTIN	08023279279030

Description

Cabinet with epoxy coated steel structure and tempered glass.
Rack-mounted adjustable glass shelves.
Rubber feet.

Made in Italy.

Technical Specifications

- Doors: 4
- Adjustable shelves: 2
- Weight: 45 kg
- Size: 75 x 38 x h 169 cm