

CABINET - 2 doors - tempered glass



Product Code	27905
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	Not required
EAN/UPC	8023279279054
GTIN	08023279279054

Description

Cabinet with epoxy coated steel structure and tempered glass.
Rack-mounted adjustable glass shelves.
Rubber feet.

Made in Italy.

Technical Specifications

- Doors: 2
- Adjustable shelves: 4
- Weight: 47 kg
- Size: 75 x 38 x h 169 cm