



Product Code	27908
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	Not required
EAN/UPC	8023279279085
GTIN	08023279279085

Description

SHELVING SYSTEM

Shelving system, made of AISI 304 stainless steel, scotch brite finish. With 5 stainless steel reinforced shelves at fixed position. Possibility to fix together several shelving systems.

Size: 100x50xh 200 cm
Load x shelf: 25 kg
Weight: 30kg

3 shelving systems fixed together

Made in Italy.