



Product Code	27911
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	Not required

EAN/UPC	8023279279115
PARAF	903608305
GTIN	08023279279115

Description

Hanging cabinet to keep medicine with 2 glass tempered doors and 2 tempered glass shelves adjustable on master rack to be hanged by two holes on the back. Supplied with lock.

Made in Italy.

Technical Specifications

- Doors: 2
- Adjustable shelves: 2
- Weight: 16 kg
- Size: 60 x 30 x h 60 cm