



Product Code	27918
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	Not required

EAN/UPC	8023279279184
GTIN	08023279279184

## Description

### GLASS-INOX CABINET

Made of AISI 304 stainless steel with two lockable doors in tempered glass 5 mm thick, 4 shelves in tempered glass 6 mm thick and 4 adjustable feet.

Size: 100x48xh 190 cm.  
Weight: 105 kg.

Made in Italy.