

PLASTIC TRAY 60x40x5h cm



Product Code	27924
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279279245

Description

Plastic tray for ISO trolleys (codes 27934-35-36).

Intelligent solutions for storing all kind of drugs, bandages, accessories and small components in a space saving and clear manner.

- Type: closed
- Dimensions: 600 x 400 x 50 mm
- Material: ABS grey