

MINI 600 ISO SERVICE TROLLEY - empty



Product Code	27934
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2552979
CND	V0899
EAN/UPC	8023279279344
GMDN	40542

Description

AURION MINI 600 - ISO trolley - empty

Size: 59 x 99.3 x h 88 cm

Basic body made in Baydur Polyurethane - Bayer AG: it can be washed and disinfected. It does not react with diluted acids, alkali and solvents. Self-extinguishing according to class UL V0.

- Latex Free.

For safety, handles are moulded into the body on both sides, trolley top has a 4-sided moulded containing rim to avoid any spilling/falling of materials.

No sharp edges.
All around bumper.
Lightweight and easy to move.

Adaptable to any model of ISO 600x400 baskets and trays thanks to the special rail sides, which allow an easy sliding and pull out.

4 wheels Ø 125 mm (3 with brake)

Made in Italy.

