



Product Code	27935
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279279351
GTIN	08023279279351

Description

ISO SERVICE TROLLEY - empty

Service trolley for "ISO" trays or baskets (not included) made of epoxy powder painted steel.
Supplied with 5 height adjustable runners for ISO trays/baskets.

Size: 65x47xh 100 cm
Weight: 12.5 kg

Pivoting castors Ø 100 mm (2 with brake)

Made in Italy.