



Product Code	27936
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279279368
GTIN	08023279279368

Description

ISO TRANSPORT TROLLEY - empty

Transport trolley for ISO trays and baskets (not included).
 Made of stainless steel AISI 304 sheet, scotch brite finishing.
 Anti-dust gasket.
 Internal side with plastic runners for ISO trays or baskets.
 Supplied with handle for transport.

External size: 52x68xh 159 cm
 Internal size: 40x60xh 129 cm

Equipped with key lock and doors kit made of transparent polycarbonate.

Antistatic castors Ø 125 mm (2 with brake)

Made in Italy.