

TRUSTEEL NON ABSORB. SUTURE gauge 6 circle 1/2
needle 44 mm - 45 cm - natural



Product Code	22794
Unit of sale	box of 48 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	2381046
CND	H0102010199
EAN/UPC	8023279227949

Description

TRUSTEEL - non absorbable monofilament steel surgical suture

Equivalent Ethicon brand: ETHISTEEL

- Gauge: 6
- Metric: 8
- Colour: Natural
- Circle: 1/2
- Needle profile: Round body blunt point
- Needle length: 44 mm
- Suture length: 45 cm

- Technical details:
- monofilament stainless steel LVM 316 grade wire
 - optimal compatibility with body tissue
 - low or no tissue reactivity
 - non-adhering and biologically inert
 - sterilized by ethylene oxide
 - shelf life of 5 years

- Applications:
- orthopedic surgeries
 - sternum closure

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.