


TRANSPARENT PLASTIC ISO DRAWER 600x400x50 mm

- closed



Product Code	27957
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279279573
---------	---------------

Description

Plastic tray for ISO trolleys (codes 27934-35-36).

Intelligent solutions for storing all kind of drugs, bandages, accessories and small components in a space saving and clear manner.

- Type: closed
- Dimensions: 600 x 400 x 50 mm
- Material: Transparent PC