

TRANSPARENT PLASTIC ISO DRAWER 600x400x100 mm - closed



Product Code	27958
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279279580
---------	---------------

Description

Plastic tray for ISO trolleys (codes 27934-35-36).

Intelligent solutions for storing all kind of drugs, bandages, accessories and small components in a space saving and clear manner.

- Type: closed
- Dimensions: 600 x 400 x 100 mm
- Material: Transparent PC