



Product Code	27960
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279279603

## Description

### SHELL MEDICINE CABINET - white

Structure in powder-coated aluminium profiles with resistant and durable melamine-faced wood particle panels.

- 8 internal adjustable shelves and lockable cabinet
- 20 unbreakable and recyclable transparent compartment trays (10 on each door).

Edged doors revolving on 270°, hinges complete with glassfilled nylon lock and handles.  
4 powder-coated steel feet at the base with adjustable plastic insoles.  
Size: 910x600xh 1980 mm