



Product Code	28006
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279280067

## Description

### LUXOR CHAIR - mechanical - blue

Chair can move up and down and rotate 360°.  
Backrest and footrest are adjustable.  
Removable armrests.

Size: 183 x 83 x h 56-76 cm  
Net weight: 72 kg  
Max load: 175 kg  
Angle adjustment of the backrest 0-80°