



Product Code 80150

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP No

CH-REP No

EAN/UPC 8023279801507

## Description

Quick release VET tourniquet for small animals.

Extra slim, ideal for cats and small dogs.

- easy to put on
- practical handling
- smooth and elastic cord
- plastic cord clamp

Multilingual instructions: GB, FR, IT, ES, PT, DE, PL, RO, GR, SE

## Technical Specifications

Cord clamp size: 43x33x8 mm

Cord size: 200 mm

Weight: 10 g