

SOEHNLE 6858 VETERINARY SCALE



Product Code	80220
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4006484005445

Description

SOEHNLE 6858 VETERINARY SCALE

Digital vet scale with large and flat stainless steel weighing platform. Designed for small and medium-sized animals. Quick and easy operation. Easy to read LCD display with 25 mm high digits and aluminium housing. Auto-off function. Hold function and taring function.

Supplied with rubber mat and multilingual manual: GB, FR, IT, DE.

Technical Specifications

Capacity: 200 kg
Sensitivity: 50 g
Display: 25 mm backlit LCD
Weighing units: kg, lb, oz, lb.oz
Size: platform 900x600xh 40 mm - indicator 215x100xh 70 mm
Power supply: DC 5V power supply (included) or 4x 1.5 V AA batteries (not included)