



Product Code	80220
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	4006484005445
---------	---------------

Description

SOEHNLE 6858 VETERINARY SCALE

Digital vet scale with large and flat stainless steel weighing platform. Designed for small and medium-sized animals. Quick and easy operation. Easy to read LCD display with 25 mm high digits and aluminium housing. Auto-off function. Hold function and taring function.

Supplied with rubber mat and multilingual manual: GB, FR, IT, DE.

Technical Specifications

Capacity: 200 kg

Sensitivity: 50 g

Display: 25 mm backlit LCD

Weighing units: kg, lb, oz, lb.oz

Size: platform 900x600xh 40 mm - indicator 215x100xh 70 mm

Power supply: DC 5V power supply (included) or 4x 1.5 V AA batteries (not included)