



| | |
|---------------|----------------|
| Product Code | 20803 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2837029 |
| CND | T0207 |
| EAN/UPC | 8023279208030 |
| PARAF | 980296800 |
| GTIN | 08023279208030 |
| GMDN | 38260 |

Description

Wide range of patterns to meet all healthcare professionals needs.
Light and comfortable, provided with adjustable laces.

Available in two sizes: medium size (M), suggested for men and women with short or medium-long hair; large size (L), suggested for women with long, curly hair.

100% high quality pure cotton, highly breathable surgical caps. MD Class I.