



Product Code	80301
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	Not required
--------	--------------

EAN/UPC	8023279803013
---------	---------------

## Description

### VETERINARY OPERATION TABLE

Veterinary operating table with electric lift (82/102 cm), baked enamelled steel tubular frame and AISI 304 stainless steel top (120x56 cm) with hole and a liquid bucket holder underneath. Manual Trendelenburg 15° and Antitrendelenburg up to 15°. Lateral hook to immobilize animals.

Made in Italy.

Size: 118x56xh 85/102

Weight: 48 kg

Max load: 90 kg