

# CLEOPATRA ELECTRIC CHAIR 3 motors - blue



Product Code	28041
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279280418

## Description

### CLEOPATRA ELECTRIC CHAIR - 3 motors - blue

Electrically controlled 3 motors chair to adjust separately height, backrest and seat inclination by foot switch or push button panel.  
Footrest adjustable by gas lift pump.  
Reversible arms and adjustable headrest.  
Cushion included.

Operating power: 100-240 V - 50/60 Hz.

Manual: GB, FR, IT, ES.

Size: 196 x 82 x h 58-100 cm  
Net weight: 88 kg  
Max. load: 175 kg  
Height adjustment: 58-100 cm  
Footrest extension: 15 cm adjustable by gas lift pump  
Face hole behind cushion